



INSTALLATION QUALIFICATION PROCEDURES

TRILUTION® micro v2.4 for Tablet



The enclosed Installation Qualification (IQ) procedures are aimed at end users who are implementing GLP-type requirements and can be readily incorporated into proprietary Standard Operating Procedures (SOPs).



1. Product Identification

Manufacturer

Gilson, Inc.
3000 Parmenter Street
PO Box 620027
Middleton, WI 53562-0027 USA

Name: TRILUTION® micro

Version: 2.4

Supplier

Organization: _____

Address: _____

Phone Number: _____ Fax Number: _____

Service Technician: _____

User

Organization: _____

Department: _____

Site (Room): _____

Date of Installation: _____

System that Includes the Software: _____

Primary Contact:

Name: _____ Phone: _____

E-Mail: _____



2. Verification

Obtain the version for each software component in TRILUTION® micro by doing the following:

1. Power the tablet on. (TRILUTION micro starts automatically on the tablet.)
2. Log in.
 - The default User name is 'admin' and the default Password is 'Gilson268'.
 - The Password is case-sensitive.
3. Select **About** in the lower right corner of the screen.
4. Select **System info**.
5. Optionally, select **Export** to save the system info to a file.

Name	Observed Version	Expected Version
Gilson TRILUTION micro		2.4.0.1
Gilson TRILUTION micro Execution Engine		1.2.41.0
Gilson TRILUTION micro Embedded Application		1.3.1.0
Gilson TRILUTION micro GBS Interface		1.4.1.0
Gilson TRILUTION micro Normalization Assistant		1.4.2.0
Gilson TRILUTION micro qPCR Assistant		1.4.2.0

Available Documents

The following documents are provided. For reference purposes, specify the storage location for each.

LT3203USB USB, TRILUTION® micro v2 Documentation

Location: _____

LT375087 TRILUTION® micro v2 Overview

Location: _____

LT255556 TRILUTION® micro v2 Validation Certificate

Location: _____

LT375090 Normalization Assistant Validation Certificate

Location: _____

LT255553 qPCR Assistant Validation Certificate

Location: _____



3. Acceptance

If no installation problems were detected or installation problems that were detected were resolved, have the person who verified the software installation provide the information requested below.

Gilson Representative: _____

Organization: _____

Signature: _____

Date: _____

End user was provided with a copy of this document.