



Identifying and Collecting Weak UV Absorbance and Low Level Peaks of Interest From a Complex Mixture Using Conditional Logic Purification on Mass and UV Signals

Application Note PHA0411

Keywords

Gilson GX-271 LC/MS Purification System, Fraction Collection, Active Pharmaceutical Ingredient (API), Metoprolol, Verapamil, Chloramphenicol, Doxepin, Sulfamethazine

Introduction

Purification of Active Pharmaceutical Ingredients (APIs) often yields high numbers of fractions, often leading to a bottleneck of fractions to dry-down and subsequently analyze. UV-based fraction collection based on slope or threshold values offers some specificity and reduction in the number of fractions collected compared to basic time or volume based collection. LC/MS-based purification provides the highest specificity as a direct result of target ion collection from the API of interest. This application note will demonstrate the improvement in specificity offered by LC/MS based purification when compared to traditional UV alone using the Gilson GX-271 LC/MS Purification System (see figure 1).



Figure 1. Gilson GX-271 LC/MS Purification System



Mass based purification offers enhanced selectivity for purifying APIs often resulting in higher purity and fewer fractions to process. There is limited information on the effectiveness of mass based purification to accomplish this goal with samples that closely resemble an API purification run. Typical synthesis steps and reactions often lead to impurities and partial products. In some cases, the compound of interest will not be the predominant compound. This represents challenges in purification with traditional UV based injection and fraction collection. Utilizing a mass signal in conjunction with conditional logic collection parameters and complimentary UV data can provide the selectivity to collect only APIs of interest. This selectivity has shown results that dramatically reduce collected fractions, collecting only the target ions of interest for an efficient purification process versus traditional UV-only based purification.

In this study, API purification was simulated using a mixture of several pharmaceuticals at varying concentration levels. The compound of interest representing the API has a low, relative percent area compared to other compounds present in the same sample solution.

Materials & Methods

Materials

Chemicals and reagents were obtained from various scientific suppliers. All solvents used were HPLC grade or higher. All reagents were ACS grade or better.

Gilson LC/MS Purification System:

- GX-271 Liquid Handler
- Direct Injection Module – 5 mL sample loop
- 322 Pumping System H2 Heads (0.25 – 30 mL/min)
- 155 UV/VIS Detector
- 307 Make Up Pump
- MRA Splitter
- Flexar SQ 300 Single Quadropole Detector
- TRILUTION LC 2.1 SP5 LCMS Software

Purification Column:

Phenomenex - Axia 21.2 mm x 50 mm; Luna 5micron C18 (2), 100 A



Mass Spectrometer (MS) & Fraction Collection Parameters:

MS Parameters:

Ionization Mode: Positive
Channel 1: Full Scan 200 – 500 amu
Channel 2: Target (EIC): m/z 267.3
Adduct 1: 1
Adduct 2: 23

Conditional Logic Fraction Collection Parameters:

Primary Conditions: UV 254 nm
Slope Collection: FS =25, BS=25, P.W. = 0.2 AND Statement 10.0 min – 95% B
Secondary Conditions: MS Target Mass 267. 3
Level: $\geq 500,000$ TIC

Methods – Sample Preparation

Sample Solution:

40 mg/mL stock solutions of four compounds (Metoprolol, Verapamil, Chloramphenicol, Doxepin, and Sulfamethazine), were mixed together in Dimethyl Sulfoxide (DMSO) resulting in ~ 8 mg/mL of each compound.






Compound	Molecular Weight	Chemical Structure
Metoprolol	267.36	
Verapamil	454.60	
Chloramphenicol	323.13	
Doxepin	279.37	
Sulfamethazine	278.33	

Table 1. Compound Solutions and Corresponding Molecular Weights



Results

Performing a fraction collection using only UV detection and peak slope collection produces seven different fractions. Each of these would be subdivided into two or three fractions resulting in ~20 total fractions consisting of 15 or 20 mL volumes. This is shown below using Fraction Simulation within TRILUTION LC to provide a 'best-fit' criteria for fraction collection.




Figure 2. Fraction Simulation Using UV and Slope Collection Mode



The simulated fraction collection for the API of interest (Metoprolol) is in relatively low amount compared to the other impurities or peaks of non-interest. The same compound mixture was purified with conditional logic collection utilizing both a UV wavelength of 254 nm and a target mass of 267.36 (Metoprolol). The Metoprolol target mass channel showed tailing. As a result, the next UV peak was eluted while the target mass of 268.2 [M+H] was still present, which caused a second fraction to be collected. The gradient conditions could be modified to increase the separation and prevent the second fraction. Increasing the make-up flow rate to 0.5 mL/min was attempted to reduce the amount of tailing with minimal improvement.




Figure 3. Gilson TRILUTION® LC Software Conditional Logic Fraction Collection Purification Based on UV and MS Signals – Fraction Tracking Function



Sample Table

Injection Number	Peak Name	Retention Time (min)	Area (uVmin x100)	Area %	Height (mV)	Sample Name	Fraction Site(s)
2	Chloramphenicol	0.767	5484120.8781	32.234	120.713	New Mix	
2	Sulfamethazine	2.988	6870992.4011	40.385	112.655	New Mix	
2	Doxepin	4.673	4310787.5011	25.337	150.053	New Mix	
2	Metroprolol	6.823	21246.6674	0.125	1.241	New Mix	Fraction Zone-1
2	Vanillin	8.367	80784.1664	0.477	3.931	New Mix	Fraction Zone-2