




# Purification of a Protein

## using the PLC 2020 Personal Purification System



### Conditions

Sample: Protein (proprietary)

Load: 0.1–10 mg

HPLC column: GE Mono Q glass column,  
10 x 100 mm

Mobile phase:

- A: 20 mM Tris buffer pH 8  
10 mM  $\beta$ -mercaptoethanol
- B: 20 mM Tris buffer pH 8  
10 mM  $\beta$ -mercaptoethanol  
1 M sodium chloride

Wavelength: 254 & 280 nm

Small scale purification of a protein was performed using anion exchange. The protein peak of interest collected was within fraction collection tube 8 (see arrow). Time per tube collection was performed to capture all UV and non-UV absorbing peaks for further research.

To discover more Application Briefs and notes from Gilson, visit [www.gilson.com/en/AI/Applications/](http://www.gilson.com/en/AI/Applications/)