

platemaster®

FAST MANUAL MICROPLATE PIPETTING



Save Time And Increase Your Productivity

With minimal setup you can fill 96- or 384- well plates in less than one minute by reducing the number of manipulations for filling a microplate.

Simplifies Sample Prep

PLATEMASTER® requires no programming, no power and minimal training, so you can start pipetting right away.

Convenient For Any Microplate Application

From plate transfers to supernatant, achieve accuracy and precision on the bench or under the hood.

PLATEMASTER is an amazingly easy-to-use, accurate solution for high throughput pipetting of 96- and 384-well microplates.

PLATEMASTER		
PLATEMASTER P20 (F110761)	0.5 μL - 20 μL Adjustable in 0.1 μL increments	
PLATEMASTER P220 (F110762)	2 μL - 220 μL Adjustable in 1 μL increments	
Dimensions	300 mm x 450 mm x 400 mm (W x L x H)	
Weight	16 kg	

PLATEMASTER ACCESSORIES		
Description	Part Number	
F1077602	384-well adapter, positioned done by wheel.	
F1077603	384-well adapter, positioned done by hand.	
F1077604	Alu-heater block for PLATEMASTER®, 96 x 0.2 mL, for PCR tubes, stripes & PCR plates Ø 7.8 mm	
F1077605	Pipetting head height adjusters (2 units)	
F1077606	Lubrication box of O-rings PLATEMASTER®	

ISO 17025 & GILSON PIPETTES

Gilson Pipette with a Calibration Certificate Issued by ISO 17025 Accredited Lab

Around the world, the accreditation of pipette testing and calibration laboratories is based on the ISO 17025 international standard. ISO 17025 accreditation enables test data to be internationally validated and accepted.

This system of international mutual recognition allows test data to be equivalent throughout the world because they are linked to the same metrological chain of international standard

weight.

Hence, the measure of the uncertainties derived from equipment, technicians, pipettes and tips which indicates the error range attributable to the measured quantities allows the accredited laboratory to attest to the pipette's reliability and ensures that it is within acceptable limits after calibration.

- Gilson ISO 17025-accredited laboratory in France provides the undisputable guarantee and confidence that our Quality System and Standard Operating Procedures have been thoroughly evaluated and approved by a legitimate third party.
- Gilson offers accredited calibration for new and serviced pipettes in full compliance with the ISO 17025 & ISO 8655 requirements.
- A certificate of calibration, approved by the French organization of accreditation (COFRAC), is delivered with each pipette. The certificate is recognized by international organizations of accreditation.
- In an on-going effort to address the growing demands for enhanced traceability and affordable ISO 17025 solutions, Gilson has been granted accreditation by the COFRAC for a special Gilson 3 volumes x 4 measurements method in addition to being accredited for the standardized ISO 8655 3 volumes x 10 measurements method. This enables us to offer ISO 17025 calibration at an attractive price.
 - The standardized ISO 8655 method (3 volumes x 10 measurements).
 - The special Gilson method (3 volumes x 4 measurements).



New Pipettes with a Calibration Certificate Issued by the Gilson ISO 17025 Accredited Lab

To order a new pipette with a calibration certificate, add a $\bf C$ at the end of the reference for 3 x 10 measurements method (standard ISO 8655), or a $\bf CG$ for a 3 x 4 measurements (Gilson method, available in French service centre recognised by accreditation organization).

Example:

- You need a P100 delivered with a 3 x 10 measurements:
 Order the following reference: F123615C
- You need a P2N delivered with a 3 x 4 measurements:
 Order the following reference: F144561CG



Please contact your nearest distributor for calibration certificate request form.





Serviced Pipettes with a Calibration Certificate Issued by ISO 17025 Accredited Lab

Some Gilson service centers in the world are accredited (ISO 17025) for pipette calibration when using the gravimetric method (ISO 8655).

If your nearest service center is accredited, you can ask for a certificate of calibration including information about the conformity of your pipette. This certificate is recognized throughout the world and linked to the International System of Units.

If your local Gilson service center is not yet accredited, Gilson's Accredited Laboratory of Calibration, based in France, can provide you the following service offers:

ISO 8655 3 X 10 MEASUREMENTS METHOD (NO TIP CHANGE)		
Part Number	Description	
PCI 001	Fixed volume Pipettes	
PCI 002	Adjustable volume Pipettes	
PCI 003	8 channel Pipettes	
PCI 004	12 channel Pipettes	

GILSON 3 X 4 MEASUREMENTS METHOD (NO TIP CHANGE)		
Part Number	Description	
PGCI 001	Fixed volume Pipettes	
PGCI 002	Adjustable volume Pipettes	
PGCI 003	8 channel Pipettes	
PGCI 004	12 channel Pipettes	